

TriDix Bio
INNOVATING DRUG DISCOVERY

PRODUCT DATA SHEET

HEK293 GCGR β -Arrestin Cell Line

SKU: TDC-S1502

Product Details

Product Description: Recombinant HEK293 cells stably overexpress human glucagon-like peptide-1 (GCGR) tagged with smBit on the surface and β -Arrestin tagged with Lgbit.

Catalog Number: TDC-S1502

Source: Genetically engineered HEK293 derivative from single clone

Package Size: Two vials of frozen cells ($>1 \times 10^6$ per vial in 1 mL)

Stability: Stable through more than 15 passages with no significant changes in assay performance or expression profile.

Culture Properties: Adherent

Associated Media: GPCR cell line culture medium kit (Cat.# TDM-1017K) containing basal media and supplements

Storage Conditions & Shipment

Product Format/Shipped: Cryopreserved / Dry ice

Storage: Liquid Nitrogen

Safety Precaution

PLEASE READ BEFORE HANDLING ANY FROZEN VIALS. *Please wear appropriate Personal Protection Equipment (lab coat, thermal gloves, safety goggles and a face shield) when handling*

Description

The HEK293 GCGR β -Arrestin Cell Line is a stable human cell line engineered to monitor β -arrestin recruitment to the glucagon receptor (GCGR). This system enables sensitive, real-time detection of ligand-induced GPCR signaling.

This stable expression system provides a robust and reproducible platform for studying GCGR signaling and for screening small molecules, peptides, or biologics targeting this clinically relevant GPCR.

Product Data




Figure 1, Dose-dependent activation of GCGR in the HEK293 GCGR β -arrestin stable cell line. Cells were stimulated with increasing concentrations of Glucagon ligand, and β -arrestin recruitment was quantified by luminescence readout. Data show a robust, ligand-dependent signal more than a 10-fold signal-to-basal ratio.

Applications

- Functional analysis of glucagon receptor signaling
- Screening of GCGR agonists, antagonists, and modulators
- Mechanistic studies of β -arrestin-mediated GCGR pathways
- Comparative pharmacological profiling in metabolic research

Key Features:

- **Engineered reporter design:** GCGR and β -arrestin are fused to complementary luciferase fragments (SmBiT/LgBiT), reconstituting active luciferase upon receptor activation and arrestin recruitment.
- **Robust and reproducible signal:** Glucagon stimulation elicits a strong, dose-dependent luminescent response with high signal-to-background ratio.
- **Assay-ready format:** Optimized for use in 96- and 384-well plate formats, enabling medium- to high-throughput screening.
- **Efficient and cost-effective:** Assay completion within 2 hours, with low material cost per plate.

- **Validated performance:** Consistent responsiveness to glucagon and related analogs across passages. Ready-to-use format: Cells are shipped as cryopreserved vials with supporting protocols for thawing, culture, and assay setup.

Protocols:

Thawing HEK293 GCGR β -Arrestin Cell

Procedures:

1. Add 30ml of GPCR cell culture supplement (Cat.# TDM-1017A) into GPCR cell line base medium (Cat.# TDM-1017).
2. Remove cells from small liquid nitrogen dewar
3. Thaw in water bath at 37°C for 30sec to 1 min.
4. Pour whole tube in 15mL Falcon tube containing 10ml of GPCR cell culture medium.
5. Spin down at 200g at room temperature for 10min.
6. Aspirate media and add 25mL medium to resuspend pellet. Pipette up and down. Make sure cells disperse well.
7. Add 12mL into each new flask (e.g. 10cm dish or T75 flask)
8. Put the vessels in 37°C 5% CO₂ incubator

Subculturing of HEK293 GCGR β -Arrestin Cell

Procedures:

1. Aspirate old media on top of cell monolayer
2. Add 5mL Ca²⁺/Mg²⁺-free Dulbecco's phosphate buffered saline (DPBS) to wash/rinse cells.
3. Aspirate DPBS and add 2mL Trypsin-EDTA solution. Leave in incubator 2-5min until cell layer is dispersed.
4. Remove plate from incubator and add 8mL supplemented cell culture medium to the flask.
5. Gently pipette up and down 5-8 times to get single cells.
6. Make a 1:5 or 1:3 dilutions and then put in 37°C, 5% CO₂ incubator

Subculturing of HEK293 GCGR β -Arrestin Cell

Procedures:

1. Culture the desired quantity of HEK293 GCGR β -Arrestin cells to 70-90% confluency
2. Remove the cells from the tissue culture flasks by washing with DPBS, adding 2mL trypsin,
3. adding 8mL media to stop the reaction, then transferring to a 15mL tube.
4. Centrifuge the cell suspension at 200 x g for 10 min at room temp.
5. Aspirate the medium and re-suspend pellet in the pre-determined volume of chilled freezing medium.
6. Determine viable and total cell counts and calculate volume of freezing medium required to

7. yield a final cell density of > 3 Million viable cells/mL
8. Dispense aliquots of this suspension into cryovials
9. Freeze cells in an automated, control-rate freezing apparatus or using a manual method
10. following standard procedures.
11. Transfer vials to liquid nitrogen storage the next morning.